

EXAMINATION GLOVE PORTFOLIO



INFECTION PREVENTION

O&M Halyard provides a full line of single-use infection control solutions to help reduce the spread of infections among patients and staff, from examination areas, to labs, to the OR[†]. Included in this range are Facial Protection Solutions, Examination Gloves, Headwear, Footwear, Gowns and Scrub Suits.




TABLE OF CONTENTS

| | |
|---|----|
| HALYARD'S NITRILE GLOVE PORTFOLIO | 4 |
| GLOVES PROPERTIES | 6 |
| CHEMOTHERAPY RESISTANCE BARRIER GUIDE | 6 |
| CHEMICAL RESISTANCE BARRIER GUIDE | 8 |
| PHYSICAL PROPERTIES | 10 |
| OTHER PRODUCT SPECIFICATIONS & QUALITY STANDARDS | 12 |
| EDUCATIONAL | 14 |
| GLOVES BY DEPARTMENT | 14 |
| SMARTPULL* DISPENSING | 16 |
| PROPER DORNING AND DOFFING OF NON-STERILE EXAM GLOVES | 17 |
| ITEM INDEX | 18 |



HALYARD'S NITRILE GLOVE PORTFOLIO

OFFERING A GLOVE FOR EVERY NEED



PURPLE NITRILE* Exam Gloves

Maximum protection for high risk tasks or situations

PURPLE NITRILE* Exam Gloves are frequently the glove choice for healthcare workers who need excellent barrier protection from bacteria, viruses and chemotherapy drugs. These durable gloves are available in 240, 310 and 400 mm lengths for maximum coverage and protection. Halyard's PURPLE NITRILE* Exam Gloves are the original PURPLE NITRILE* Gloves produced by Safeskin.



STERLING* Nitrile Exam Gloves

Moderate risk tasks require a standard glove with STERLING* reliability

STERLING* Nitrile exam gloves were specifically designed to fit and feel like latex, with excellent comfort and tactile sensitivity. These gloves provide durable barrier protection and can be used facility-wide in a multitude of tasks. These gloves are available in 240 and 310 mm lengths for greater coverage and protection.



BASICS* Blue & BASICS* White Nitrile Exam Gloves

Low to moderate risk tasks require an economical glove that does not compromise on protection or comfort

BASICS* Blue & BASICS* White Nitrile exam gloves are thinner, lighter, and more economical than our standard nitrile gloves. The gloves are suitable for a wide range of tasks and are the ideal replacement for vinyl or latex exam gloves while maintaining comfort and fitting barrier protection during use.



SKYBREEZE* Nitrile Exam Gloves

Designed to be comfortable, economical, and protective

SKYBREEZE* Nitrile exam gloves are a versatile option for various tasks and serve as a great alternative to latex or vinyl exam gloves, providing both comfort and reliable barrier protection. The gloves are approved for use with numerous chemotherapy drugs, for better peace of mind when coming into contact with these potentially hazardous drugs.



BLACK-FIRE* Nitrile Exam Gloves

The functionality of two gloves in one

BLACK-FIRE* Nitrile Exam gloves are reversible and can be used as a black glove or turned inside out to reveal the high-visibility, lower friction orange side.

INCREASING PROTECTION

PURPLE NITRILE* Outstanding barrier protection

Latex-Free
Powder-Free
Non-Sterile/Sterile
Heavy-Weight

Barrier Protection
and enhanced
performance

>240 mm/
310 mm/400 mm
length



STERLING* Nitrile Facility-wide use

Latex-Free
Powder-Free
Non-Sterile
Medium-Weight

Protection with
enhanced tactile
sensitivity

>240 mm / 310 mm
length



SKYBREEZE* Light and comfortable

Latex-Free
Powder-Free
Non-Sterile
Light-Weight

Barrier Protection
Comfortable
Cost-effective

>240 mm length



BASICS* Blue & BASICS* White Nitrile protection and comfort

Latex-Free
Powder-Free
Non-Sterile
Light-Weight

Economical
Comfortable
Protective

>240 mm length



GLOVES PROPERTIES

CHEMOTHERAPY RESISTANCE BARRIER GUIDE



The HALYARD* range of Nitrile Exam Gloves are designed to offer protection while working with chemotherapy drugs. The gloves are available in three lengths, providing protection and coverage up to the mid-forearm. The **PURPLE NITRILE* Examination glove** family are extensively tested for use with chemotherapy drugs. Our range of PURPLE NITRILE* Examination Gloves have been tested according to the European norm EN ISO 374-1 and the US standard ASTM D6978 to ensure your safety when coming into contact with potentially hazardous chemotherapy drugs.

| Chemotherapy Drugs | PURPLE NITRILE* | PURPLE NITRILE-XTRA* | PURPLE NITRILE MAX* |
|-------------------------------|----------------------|----------------------|---------------------|
| | Glove length >240 mm | Glove length 310 mm | Glove length >400mm |
| Chemotherapy Drugs | Breakthrough rating | Breakthrough rating | Breakthrough rating |
| Arsenic Trioxide (1 mg/ml) | ✓ | ✓ | ✓ |
| Azacitidine (25 mg/ml) | ✓ | ✓ | ✓ |
| Bendamustine (5 mg/ml) | ✓ | ✓ | ✓ |
| Bortezomib (1 mg/ml) | ✓ | ✓ | ✓ |
| Bleomycin sulfate (15 mg/ml) | ✓ | ✓ | ✓ |
| Busulfan (6 mg/ml) | ✓ | ✓ | ✓ |
| Carboplatin (10 mg/ml) | ✓ | ✓ | ✓ |
| Carfilzomib (10 mg/ml) | ✓ | ✓ | ✓ |
| Carmustine (3.3 mg/ml) | ✗ | ✓ | ✓ |
| Cetuximab (Erbixux) (2 mg/ml) | ✓ | ✓ | ✓ |
| Cisplatin (1 mg/ml) | ✓ | ✓ | ✓ |
| Cyclophosphamide (20 mg/ml) | ✓ | ✓ | ✓ |
| Cytarabine HCL (100 mg/ml) | ✓ | ✓ | ✓ |
| Cytovene (10 mg/ml) | NOT TESTED | ✓ | ✓ |
| Dacarbazine (10 mg/ml) | ✓ | ✓ | ✓ |
| Daunorubicin HCL (5 mg/ml) | ✓ | ✓ | ✓ |
| Decitabine (5 mg/ml) | ✓ | ✓ | ✓ |
| Docetaxel (10 mg/ml) | ✓ | ✓ | ✓ |
| Doxorubicin HCL (2 mg/ml) | ✓ | ✓ | ✓ |
| Ellence (2 mg/ml) | ✓ | ✓ | ✓ |
| Eribulin Mesylate (0.5 mg/ml) | ✓ | ✓ | ✓ |
| Etoposide (20 mg/ml) | ✓ | ✓ | ✓ |
| Fludarabine (25 mg/ml) | ✓ | ✓ | ✓ |
| Fluorouracil (50 mg/ml) | ✓ | ✓ | ✓ |
| Fulvestrant (50 mg/ml) | ✓ | ✓ | ✓ |
| Gemcitabine (38 mg/ml) | ✓ | ✓ | ✓ |
| Idarubicin (1 mg/ml) | ✓ | ✓ | ✓ |
| Ifosfamide (50 mg/ml) | ✓ | ✓ | ✓ |
| Irinotecan (20 mg/ml) | ✓ | ✓ | ✓ |
| Mechlorethamine HCL (1 mg/ml) | ✓ | ✓ | ✓ |
| Melphalan (5 mg/ml) | ✓ | ✓ | ✓ |



| | PURPLE NITRILE* | PURPLE NITRILE-XTRA* | PURPLE NITRILE MAX* |
|-------------------------------|----------------------|----------------------|---------------------|
| | Glove length >240 mm | Glove length 310 mm | Glove length >400mm |
| Chemotherapy Drugs | Breakthrough rating | Breakthrough rating | Breakthrough rating |
| Methotrexate (25 mg/ml) | ☑ | ☑ | ☑ |
| Mitomycin (0.5 mg/ml) | ☑ | ☑ | ☑ |
| Mitoxantrone (2 mg/ml) | ☑ | ☑ | ☑ |
| Oxaliplatin (2 mg/ml) | ☑ | ☑ | ☑ |
| Paclitaxel (6 mg/ml) | ☑ | ☑ | ☑ |
| Paraplatin (10 mg/ml) | ☑ | ☑ | ☑ |
| Pemetrexed (25 mg/ml) | ☑ | ☑ | ☑ |
| Pertuzumab (30 mg/ml) | ☑ | ☑ | ☑ |
| Raltitrexed (0.5 mg/ml) | NOT TESTED | ☑ | ☑ |
| Retrovir (10 mg/ml) | NOT TESTED | ☑ | ☑ |
| Rituximab (10 mg/ml) | ☑ | ☑ | ☑ |
| Temsirolimus (25 mg/ml) | ☑ | ☑ | ☑ |
| ThioTEPA (10 mg/ml) | ☒ | ☑ | ☑ |
| Topotecan HCL (1 mg/ml) | ☑ | ☑ | ☑ |
| Trastuzumab (21 mg/ml) | ☑ | ☑ | ☑ |
| Triclosan (1 mg/ml) | NOT TESTED | ☑ | ☑ |
| Trisenox (1 mg/ml) | ☑ | NOT TESTED | ☑ |
| Vinblastine (1 mg/ml) | ☑ | ☑ | ☑ |
| Vincristine Sulfate (1 mg/ml) | ☑ | ☑ | ☑ |
| Vinorelbine (10 mg/ml) | ☑ | ☑ | ☑ |
| Zoledronic Acid (0.8 mg/ml) | ☑ | ☑ | ☑ |

Rating Code

| | |
|-------------------|--|
| ☒ NOT RECOMMENDED | Avoid using this glove with this chemical |
| ⚠ CAUTION | This glove is suitable for this chemical under careful control of its use (breakthrough time < 59 min) |
| ☑ RECOMMENDED | This glove is well suited for use with this chemical (breakthrough time > 60 min) |
| NOT TESTED | - |

PURPLE NITRILE* : The following chemotherapy drugs and concentrations showed breakthrough detected in less than 30 minutes. Carmustine (3.3 mg/ml): 1.8 minutes ThioTEPA (10.0 mg/ml): 1.6 minutes. Warning: Not for use with Carmustine and ThioTEPA
 Caution: The testing conditions used are intended to approximate the worst case conditions for clinical use. Testing was conducted on single layer glove material. It is the users' responsibility to adapt the applicability of the these gloves for their intended use with the cited chemotherapy drugs. Breakthrough time data available on request.

CHEMICAL RESISTANCE BARRIER GUIDE



The **HALYARD* STERLING* Nitrile** and **PURPLE NITRILE* Examination Gloves** are cleared for use with various chemical substances according to EN ISO 374-1 (EN 16523-1) and ASTM F739-12, thereby ensuring healthcare workers are always informed and adequately protected when working with potentially hazardous chemicals.

| | STERLING* Nitrile Glove length >240 mm | PURPLE NITRILE* Glove length >240 mm | PURPLE NITRILE-XTRA* Glove length 310 mm |
|---|---|---|---|
| Chemical | Breakthrough rating | Breakthrough rating | Breakthrough rating |
| Adhesive Tape Remover Pads | ⚠ | ✔ | ✔ |
| Alcoholic Bouin's Fixative | ⚠ | ✔ | ✔ |
| Alcoholic Eosin | ⚠ | ⚠ | ⚠ |
| Bleach | ✔ | ✔ | ✔ |
| Bouin's Fixative | ✔ | ✔ | ✔ |
| B-Plus Fixative | ✔ | ✔ | ✔ |
| CarNoy's Fluid | ⊘ | ⊘ | ⊘ |
| Chlorhexidine Gluconate 4% | NOT TESTED | ✔ | ✔ |
| Cidex 14-Day | ✔ | ✔ | ✔ |
| Cidex OPA | ✔ | ✔ | ✔ |
| Decalcifier I | ✔ | ✔ | ✔ |
| EM Fixative | ✔ | ✔ | ✔ |
| Eosin Stain | ⚠ | ✔ | ✔ |
| Ethidium Bromide (1%) | ✔ | ✔ | ✔ |
| Ethyl Alcohol (95%) | ⚠ | ⚠ | ⚠ |
| Formalin Alcohol (tested for Formaldehyde) | ⚠ | ⚠ | ✔ |
| Formalin Alcohol (tested for Reagent Alcohol) | ⚠ | ✔ | ✔ |
| Formalin, 10% Buffered (tested for Formaldehyde) | ✔ | ✔ | ✔ |
| Formalin, 10% Buffered (tested for Methanol) | ✔ | ✔ | ✔ |
| Giemsa Stain Solution | ⚠ | ⚠ | ⚠ |
| Glutaraldehyde 4% | ✔ | ✔ | ✔ |
| Haematoxylin Stain (Harris Solution) | ✔ | ⚠ | ✔ |
| Hydrochloric Acid, 35.5% | ✔ | ✔ | ✔ |
| Hydrogen Peroxide, 3% | ✔ | ✔ | ✔ |
| Isopropyl Alcohol, 70% | ⚠ | ⚠ | ✔ |

| | STERLING* Nitrile | PURPLE NITRILE* | PURPLE NITRILE-XTRA* |
|---|----------------------|----------------------|----------------------|
| | Glove length >240 mm | Glove length >240 mm | Glove length 310 mm |
| Chemical | Breakthrough rating | Breakthrough rating | Breakthrough rating |
| Methanol | ⚠ | ⚠ | ⚠ |
| Methyl Methacrylate (Bone Cement) | ⊘ | ⊘ | ⚠ |
| Monsef's Solution | ✓ | ✓ | ✓ |
| OxyCide Concentrate | ⚠ | ✓ | ✓ |
| OxyCide Ready to Use | ✓ | ✓ | ✓ |
| Permaslip Mounting Medium & Liquid Cover Slip | ⊘ | ⊘ | ⊘ |
| Rapid Bone Decalcifier | ✓ | ✓ | ✓ |
| Reagent Alcohol | ⚠ | ⚠ | ⚠ |
| Resert XL HLD | ⚠ | ✓ | ✓ |
| Sodium Hydroxide (40%) | ✓ | ✓ | ✓ |
| Trichloroacetic Acid, (10%) | ✓ | ✓ | ✓ |
| Xylene, 99% | ⊘ | ⊘ | ⊘ |
| Acetone | ⚠ | ⚠ | ⚠ |
| Ethyl Acetate | ⚠ | ⚠ | ⚠ |
| n-Heptane | ✓ | ✓ | ✓ |
| Sulphuric Acid (96%) | ⚠ | ⚠ | ⚠ |
| Nitric acid (65%) | ⚠ | ⚠ | ⚠ |
| Acetic acid (99%) | ⚠ | ⚠ | ⚠ |
| Ammonium Hydroxide (25%) | ⚠ | ⚠ | ⚠ |
| Hydrogen Peroxide (30%) | ⚠ | ✓ | ✓ |
| Formaldehyde (37%) | ✓ | ✓ | ✓ |
| Incidin Rapid (Isopropanol) | NOT TESTED | ✓ | NOT TESTED |
| Incidin Foam (Isopropanol & Ethanol) | NOT TESTED | ✓ | NOT TESTED |
| Incidin Plus (Diethylene Glycol Butyl Ether & phenoxyethanol) | NOT TESTED | ✓ | NOT TESTED |


Rating Code

| | |
|-------------------|--|
| ⊘ NOT RECOMMENDED | Avoid using this glove with this chemical |
| ⚠ CAUTION | This glove is suitable for this chemical under careful control of its use (breakthrough time < 59 min) |
| ✓ RECOMMENDED | This glove is well suited for use with this chemical (breakthrough time > 60 min) |
| NOT TESTED | - |

These gloves are not intended for applications involving direct exposure to harsh chemicals, where heavy-duty industrial gloves are required. It is the users' responsibility to adapt the applicability of these gloves for their intended use with the cited chemicals. Breakthrough time data available on request.


PHYSICAL PROPERTIES

Force at break is an indication of glove strength. It measures the force (Newton) required in order to break a glove by stretching it. EN 455-2:2015 requires that nitrile examination gloves meet a minimum force at break of 6N. The HALYARD* STERLING* Nitrile & PURPLE NITRILE* Examination Glove families all exceed the minimum force at break requirement. Our BASICS* Blue, BASICS* White and SKYBREEZE* Nitrile Examination Gloves comply to the norm with a 6N force at break.




AQL (Acceptable Quality Level) is an indicator of glove quality. The lower the AQL, the higher the product quality and lower the incidence of pinholes and other defects, which can create a barrier breach and ultimately lead to infection.

EN 455-1:2020 describes a method in which randomly sampled gloves are subjected to a water tightness test, i.e. the ability of a glove to hold 1000ml of water without leaking. For compliance to the EN 455-1:2020, the Acceptance Quality Limit (AQL) should be 1.5. The HALYARD* STERLING* Nitrile & PURPLE NITRILE* Examination Glove families all exceed (in terms of performance), the Acceptance Quality Limit with an AQL 1.0. Our BASICS* Blue, BASICS* White and SKYBREEZE* Nitrile Examination Gloves maintain an AQL 1.5.









Palm Thickness is an attribute that can be used in determining both tactile sensitivity and barrier protection. Consistency for this property is key for both durability and chemical permeation protection. The below chart indicates the palm thickness of HALYARD® gloves per glove family.





OTHER PRODUCT SPECIFICATIONS & QUALITY STANDARDS

Gauge thickness measurements

| |  |  |  |  |
|----------------|---|---|---|---|
| | PURPLE NITRILE* | PURPLE NITRILE-XTRA* | PURPLE NITRILE* MAX | STERLING* Nitrile |
| Middle finger | 0.15 mm | 0.15 mm | 0.28 mm | 0.09 mm |
| Palm | 0.13 mm | 0.12 mm | 0.24 mm | 0.08 mm |
| Cuff | 0.09 mm | 0.11 mm | 0.15 mm | 0.06 mm |
| Average length | > 240 mm | 310 mm | 415 mm | > 240 mm |

Product specifications - Palm width

| |  |  |  |  |
|---------|---|---|---|---|
| | PURPLE NITRILE* | PURPLE NITRILE-XTRA* | PURPLE NITRILE* MAX | STERLING* Nitrile |
| X-Small | 70 mm | 70 mm | - | 70 mm |
| Small | 80 mm | 80 mm | 80 mm | 80 mm |
| Medium | 95 mm | 95 mm | 100 mm | 95 mm |
| Large | 110 mm | 110 mm | 111 mm | 110 mm |
| X-Large | 120 mm | 120 mm | 121 mm | 120 mm |

Quality standards

| | EN 455-1:2020 +A1:2020 EN 455-2:2015 EN 455-3:2015 EN 455-4:2009 | EN ISO 374-1:2016 Type B | EN ISO 21420:2020 | Product CE marked as per (EU) 2016/425 regulation on Personal Protective Equipment. | Approved for food handling | Product CE marked as per Regulation (EU) 2017/745 on Medical Devices | Quality system (manufacturing sites) |
|-------------------------|---|--------------------------|-------------------|---|----------------------------|--|--------------------------------------|
| PURPLE NITRILE* | ✓ | ✓ | ✓ | Class of the equipment: PPE Category III | - | Class of device: I | ISO 13485 |
| PURPLE NITRILE-XTRA* | ✓ | ✓ | ✓ | Class of the equipment: PPE Category III | - | Class of device: I | ISO 13485 |
| PURPLE NITRILE* MAX | ✓ | ✓ | ✓ | Class of the equipment: PPE Category III | - | Class of device: I | ISO 13485 |
| PURPLE NITRILE* Sterile | ✓ | - | - | - | - | Class of device: I - sterile | ISO 13485 |
| STERLING* Nitrile | ✓ | ✓ | ✓ | Class of the equipment: PPE Category III | ✓ | Class of device: I | ISO 13485 |

Gauge thickness measurements



STERLING®
Nitrile-XTRA



BASICS® Blue &
BASICS® White



BLACK-FIRE®

| | STERLING® Nitrile-XTRA | BASICS® Blue & BASICS® White | BLACK-FIRE® |
|----------------|---------------------------|---------------------------------|-------------|
| Middle finger | 0.11 mm | 0.07 mm | 0.14 mm |
| Palm | 0.08 mm | 0.05 mm | 0.10 mm |
| Cuff | 0.07 mm | 0.04 mm | 0.09 mm |
| Average length | 310 mm | > 240 mm | > 240 mm |

Product specifications - Palm width



STERLING
Nitrile-XTRA®



BASICS® Blue &
BASICS® White



BLACK-FIRE®

| | STERLING Nitrile-XTRA® | BASICS® Blue & BASICS® White | BLACK-FIRE® |
|---------|---------------------------|---------------------------------|-------------|
| X-Small | 70 mm | 75 mm | 70 mm |
| Small | 80 mm | 85 mm | 80 mm |
| Medium | 95 mm | 95 mm | 95 mm |
| Large | 110 mm | 106 mm | 110 mm |
| X-Large | 120 mm | 116 mm | 120 mm |

Quality standards

| | EN 455-1:2020 +A1:2020 EN 455-2:2015 EN 455-3:2015 EN 455-4:2009 | EN ISO 374-1:2016 Type B | EN ISO 21420:2020 | Product CE marked as per (EU) 2016/425 regulation on Personal Protective Equipment. | Approved for food handling | Product CE marked as per Regulation (EU) 2017/745 on Medical Devices | Quality system (manufacturing sites) |
|---------------------------------|--|--------------------------------|----------------------|--|----------------------------------|---|--|
| STERLING Nitrile-XTRA® | ✓ | ✓ | ✓ | Class of the equipment: PPE Category III | ✓ | Class of device: I | ISO 13485 |
| BASICS® Blue & BASICS® White | ✓ | - | ✓ | Class of the equipment: PPE Category III | ✓ | Class of device: I | ISO 13485 |
| BLACK-FIRE® | ✓ | ✓ | ✓ | Class of the equipment: PPE Category III | ✓ | Class of device: I | ISO 13485 |

EDUCATIONAL

GLOVES BY DEPARTMENT

Choosing the correct glove based on the risk of the anticipated task is crucial in a healthcare worker's daily activities.

Chemotherapy administration - Critical care - Ambulance

High risk tasks typically require a high risk glove

| SKYBREEZE* | BASICS* Blue & BASICS* White | STERLING* Nitrile | STERLING* Nitrile-XTRA* | PURPLE NITRILE* | PURPLE NITRILE-XTRA* | PURPLE NITRILE* Sterile |
|------------|------------------------------|-------------------|-------------------------|-----------------|----------------------|-------------------------|
| ✓ | - | ✓ | ✓ | ✓ | ✓ | - |

Labour/OBGYN - Dermatology - Wound care - Minor non-invasive procedures - Burn units

High risk tasks typically require a high risk glove

| SKYBREEZE* | BASICS* Blue & BASICS* White | STERLING* Nitrile | STERLING* Nitrile-XTRA* | PURPLE NITRILE* | PURPLE NITRILE-XTRA* | PURPLE NITRILE* Sterile |
|------------|------------------------------|-------------------|-------------------------|-----------------|----------------------|-------------------------|
| - | - | - | - | - | - | ✓ |

Facility-wide use - ICU - OR - Pharmacies - Laboratories

Moderate risk tasks may require a standard glove with Sterling* reliability

| SKYBREEZE* | BASICS* Blue & BASICS* White | STERLING* Nitrile | STERLING* Nitrile-XTRA* | PURPLE NITRILE* | PURPLE NITRILE-XTRA* | PURPLE NITRILE* Sterile |
|------------|------------------------------|-------------------|-------------------------|-----------------|----------------------|-------------------------|
| ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | - |

Phlebotomy - General patient care & transportation - Routine oral care

Low risk tasks often require a cost effective glove that does not compromise on protection or comfort

| SKYBREEZE* | BASICS* Blue & BASICS* White | STERLING* Nitrile | STERLING* Nitrile-XTRA* | PURPLE NITRILE* | PURPLE NITRILE-XTRA* | PURPLE NITRILE* Sterile |
|------------|------------------------------|-------------------|-------------------------|-----------------|----------------------|-------------------------|
| ✓ | ✓ | ✓ | ✓ | ✓ | - | ✓ |

The tasks listed serve as examples of typical use, however, it is ultimately up to the end user to identify and choose the most appropriate glove for the task they will undertake.



SMARTPULL* DISPENSING



The HALYARD* SMARTPULL* system features two openings to allow for better glove dispensing and less waste!



38.2% WASTE
REDUCTION

The system features two separate openings on the box. The first, smaller opening is used when the box is full to reduce multiple dispensing, which often results in exam glove waste when gloves become unusable after falling on the floor. When the box is half empty, the second, larger opening allows for easier access to the gloves.

The patented dispensing technology is featured on our PURPLE NITRILE*, STERLING* Nitrile, BASICS* Blue & BASICS* White Nitrile and SKYBREEZE* Examination Glove boxes.

DIRECTIONS FOR USE




STEP 1

Remove the perforated carton opening at Tab 1



The smaller opening helps prevent multiple dispensing and glove waste.



STEP 2

As the box empties and it becomes difficult to reach inside, remove the opening indicated by Tab 2



Continue dispensing as you would from a traditional glove carton.

The SMARTPULL* dispensing system is available on the standard (240mm) length glove boxes.

PROPER DONNING AND REMOVAL OF NON-STERILE EXAM GLOVES

DONNING



01 Before donning your gloves, always wash your hands.



02 Take out a glove from its original dispenser, holding only the cuff.



03 Hold the glove at its opening with one hand and slide your fingers and thumb of the opposite hand into the glove.



04 Pull the glove towards your wrist to fully don the glove. Only use your knuckles while doing this to avoid fingernail punctures.



05 Take a second glove out of the dispenser with your bare hand, holding only the cuff.



06 Hold the glove at its opening and slide your fingers and thumb into the glove. Pull the glove towards your wrist using the knuckles of your gloved hand.

REMOVAL



01 Using your dominant hand. Grasp the outside cuff of your opposite gloved hand and peel the glove off. This turns the glove inside out.



02 Hold the removed glove in your gloved hand.



03 Slide two fingers of your ungloved hand under the remaining glove cuff at your wrist.



04 Peel the second glove off your hand turning it inside out and enclosing the first glove within it.



05 Discard your gloves in the appropriate waste container or garbage can.



06 After you dispose of your gloves, wash your hands thoroughly.

ITEM INDEX

| ITEM INDEX | | | | | | |
|--|---------------------------------|-------------------|--------|--------|--------|---------|
| Description | Packaging | Item code by size | | | | |
| | | X-Small | Small | Medium | Large | X-Large |
| PURPLE NITRILE* Exam Gloves | 100 gloves/box 10 boxes/case | 52000M | 52001M | 52002M | 52003M | -- |
| | 90 gloves/box 10 boxes/case | -- | -- | -- | -- | 52004M |
| PURPLE NITRILE-XTRA* Exam Gloves | 50 gloves/box 10 boxes/case | 50600M | 50601M | 50602M | 50603M | 50604M |
| PURPLE NITRILE* MAX | 50 gloves/box 8 boxes/case | -- | 47195 | 47196 | 47197 | 47198 |
| PURPLE NITRILE* Sterile Exam Gloves | 50 gloves/box 4 boxes/case | -- | 52201M | 52202M | 52203M | -- |
| STERLING* Nitrile Exam Gloves | 200 gloves/box 10 boxes/case | 13938 | 13940 | 13941 | 13942 | -- |
| | 170 gloves/box 10 boxes/case | -- | -- | -- | -- | 13943 |
| STERLING* Nitrile-XTRA* Exam Gloves | 100 gloves/box 10 boxes/case | 44286 | 44287 | 44288 | 44289 | 44290 |
| BASICS* Blue Nitrile Exam Gloves | 200 gloves/box 10 boxes/case | 44749 | 44750 | 44751 | 44752 | -- |
| | 170 gloves/box 10 boxes/case | -- | -- | -- | -- | 44753 |
| BASICS* White Nitrile Exam Gloves | 200 gloves/box 10 boxes/case | 44740 | 44741 | 44742 | 44743 | -- |
| | 170 gloves/box 10 boxes/case | -- | -- | -- | -- | 44744 |
| SKYBREEZE* | 200 gloves/box 10 boxes/case | 47373 | 47374 | 47375 | 47376 | -- |
| | 180 gloves/box 10 boxes/case | -- | -- | -- | -- | 47377 |
| BLACK-FIRE* | 150 gloves/box 10 boxes/case | 45256 | 45257 | 45258 | 45259 | -- |
| | 140 gloves/box 10 boxes/case | -- | -- | -- | -- | 45260 |



† Operating Room

For more information,
please contact your Halyard Customer Care Team:

UK & IE (toll-free)

Tel: 0800 085 82 04 - Fax: 0800 066 91 06
customerservice.uk.ie@hyh.com

South Africa (toll-free)

Tel: 0801 111 332 - Fax: 0801 111 335
customerservice.za@hyh.com

Other countries

Tel: +32 2 585 51 14 - Fax: +32 2 897 00 95
customerservice.export@hyh.com

* Registered Trademark or Trademark of O&M Halyard or its affiliates.
The colour GREY is a Registered Trademark or Trademark of O&M Halyard or its affiliates. The colour PURPLE is a Registered Trademark or Trademark of O&M Halyard or its affiliates. The colour ORANGE is a Registered Trademark or Trademark of O&M Halyard or its affiliates.
©2018, 2023. All rights reserved. HC288-04-UK COPY-01741