

EGIS

EGIS

Colorectal
Stent

EGIS **Colorectal** Stent

Large Intestine or Rectum Stent



Double Covered Type

EGIS

Colorectal Stent


EGIS Colorectal Stent

Large Intestine or Rectum Stent

Colorectal Stent



EGIS Colorectal Stent is used to open strictures in the large intestine. The Stent is designed with outstanding conformability to fit naturally into the large intestine. Colorectal Stent fits naturally into the lower digestive tract and allows patients to defecate normally. The proximal, middle and distal GOLD markers increase visibility in fluoroscopy.



Specification

Delivery System **10Fr, 80/150/220cm** | **12Fr, 80/150cm**
 Ordering No. First two** : Diameter(mm) / Last two** : Length(cm)

TYPE	Ordering No.		Size	
	OTW	TTS	Diameter	Length
Single Bare	CSB****08	CSB****22	18, 20, 22, 24, 26, 28, 30	6, 8, 10, 12
	CSB****15			
Double Bare	CDB****08	CDB****22		
	CDB****15			
Single Covered	CSC****08	CSC****22		
	CSC****15			
Double Covered	CDC****08	CDC****22		
	CDC****15			