


Trap it

TRAPPING BALLOON

...CONVENIENT, SIMPLE, SAFE




IMDS

INTERVENTIONAL MEDICAL
DEVICE SOLUTIONS

Trapit

TRAPPING BALLOON



* Only the TRP10015 has a depth marker for a 90 cm guide catheter length


BALLOON DESIGNED TO TRAP

HIGHLY PUSHABLE SHAFT

UNIQUE AUTO STOP



NOMINAL PRESSURE **8 ATM**



RATED BURST PRESSURE **12 ATM**

Model	Balloon OD at nominal	Balloon length	Shaft OD	Guiding catheter compatibility
TRP9015 (T-90)	2.5 mm	15 mm	2.0F (0.66 mm)	90 cm
TRP10015 (T-100)	2.5 mm	15 mm	2.0F (0.66 mm)	100 cm

for more information
WWW.IMDS.NL

CE
0088

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.
Caution: Federal Law (U.S.) restricts this device to sale by or on the order of a physician.
TrapIT™ is a trademark of IMDS.

Ceintuurbaan Noord 150
9301 NZ Roden
The Netherlands
sales@imds.nl
Phone +31(0)50 8200230
Fax +31(0)50 8200231

