



**A New Design**  
Reducing the Risk of AV Block

# 1 Reducing the Risk of AV Block

Elongatable waist designed for a variety of VSD

Soft woven mesh provides high conformability to septal defects

# 2 Low Profile

Partnered with 5-7 Fr. delivery sheath


Slim cable to minimize unwanted damage

# 3 Procedural Flexibility

Double-sided screw for retrograde or antegrade approach

## Product Specification

Specification		D Disc Diameter (mm)	D2 Waist Diameter LV Side (mm)	D1 Waist Diameter RV Side (mm)	L Waist Length (mm)	Recommended Delivery Sheath (Fr.)
	LT-MFO-5-3	10	5	3	4	SFP5F
	LT-MFO-6-4	10	6	4	4	
	LT-MFO-7-5	12	7	5	4	
	LT-MFO-8-6	12	8	6	4	
	LT-MFO-9-7	14	9	7	4	SFP6F
	LT-MFO-10-8	14	10	8	4	
	LT-MFO-12-10	16	12	10	4	SFP7F
	LT-MFO-14-12	18	14	12	4	



## LifeTech Scientific Corporation

Add: Cybio Electronic Building, Langshan 2nd Street, North Area of High-tech Park, Nanshan District, Shenzhen 518057, P. R. China

Tel: +86-755-86026250 Fax: +86-755-86026251

Email: marketing@lifetechmed.com

Web: www.lifetechmed.com