

MEET YOUR STERILISATION POUCH NEEDS WITH CONFIDENCE

The HALYARD* portfolio offers a line of high-quality single-use sterilisation pouches and reels to meet your sterile processing department's sterilisation packaging needs.

- Feature triple heat-sealed rail seals
- CE marked and comply with the international norms and standards – ISO 11607 & EN 868
- Allow for aseptic presentation of contents

- Available in:
 - A wide range of sizes
 - Pre-Cut Heat-Seal Pouches
 - Heat Seal Reels (Rolls)
- Feature:
 - Pre-Cut pouches feature heat-sealed corners on the thumb notch end, which help prevent dust from collecting at the opening seal and reducing the risk of contamination when pouch is opened
 - Thumb Notch for easy opening
 - Specialty designed peel-off seal that makes opening easy, clean, & fiber-free



HEAT-SEAL POUCHES AND REELS

HALYARD* Paper-Film Sterilisation Pouches are made from premium 68 GSM paper consisting entirely of primary virgin fiber with excellent bacterial barrier properties that allow for proper sterilant penetration. The product provides clean peeling performance and high seal strength. Multilayer transparent film provides puncture and abrasion resistance while also delivering a high-quality seal. The combination of our paper and film create a high-performance, durable sterile barrier system that provides users with a reliable clean peel experience. HALYARD* Paper/Film Pouches are printed with sterilisation process indicators for steam and EO sterilisation modalities.



HALYARD* Heat-Seal Paper/Film Pouches

Product Code	Size	Compatible Sterilisation Modality	Packaging
48516	50 x 250mm (1.97 x 9.84 in)	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48518	75 x 200mm (2.95 x 7.87 in)	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48495	100 x 150mm (3.94 x 5.91 in)	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48498	100 x 270mm (3.94 x 10.63 in)	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48499	100 x 300mm (3.94 x 11.81 in)	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48500	100 x 380mm (3.94 x 14.96 in)	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48501	100 x 400mm (3.94 x 15.75 in)	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48503	120 x 400mm (4.72 x 15.75 in)	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48506	150 x 300mm (5.91 x 11.81 in)	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48558	150 x 350mm (5.91 x 13.78 in)	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48508	150 x 400mm (5.91 x 15.75 in)	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48510	205 x 400mm (8.07 x 15.75 in)	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48511	220 x 300mm (8.66 x 11.81 in)	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48512	250 x 380mm (9.84 x 14.96 in)	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48513	250 x 500mm (9.84 x 19.69 in)	STEAM/EO	250 eaches/box 2 boxes/case 500 eaches/case
48514	300 x 500mm (11.81 x 19.69 in)	STEAM/EO	250 eaches/box 2 boxes/case 500 eaches/case
48515	420 x 500mm (16.54 x 19.69 in)	STEAM/EO	250 eaches/box 2 boxes/case 500 eaches/case

HALYARD* Heat-Seal Paper/Film Gusseted Pouches/Reels

48492	150 x 50 x 400m (5.91 x 1.97 x 15.75 in)	STEAM/EO	100 eaches/box 10 boxes/case 1,000 eaches/case
48494	150mm x 50mm x 100m (5.91 in x 1.96 in x 328 ft)	STEAM/EO	Shrink-wrapped individually 1 reels/case

HALYARD* Heat-Seal Paper/Film Reels

Product Code	Size	Compatible Sterilisation Modality	Packaging
48526	75mm x 200m (2.95 in x 656.2 ft)	STEAM/EO	Shrink-wrapped individually 2 reels/case
48521	100mm x 200m (3.94 in x 656.2 ft)	STEAM/EO	Shrink-wrapped individually 1 reels/case
48522	150mm x 200m (5.91 in x 656.2 ft)	STEAM/EO	Shrink-wrapped individually 1 reels/case
48523	200mm x 200m (7.87 in x 656.2 ft)	STEAM/EO	Shrink-wrapped individually 1 reels/case
48524	250mm x 200m (9.84 in x 656.2 ft)	STEAM/EO	Shrink-wrapped individually 1 reels/case
48555	300mm x 200m (11.81 in x 656.2 ft)	STEAM/EO	Shrink-wrapped individually 1 reels/case



Pouch Composition: Paper Backer - Green Tint Film Front