

UnoMeter™ Exchange Bags



Intended use

UnoMeter™ Exchange Bags are intended for use with the UnoMeter™ 500, a mechanical and disposable urine meter which is intended to be used for collection, sampling and measuring urine volume output (diuresis) over time in critically ill patients. The products are intended for single use only.

Exchange bag 2L:



Exchange bag 3L:



UnoMeter™ Exchange Bags:

REF No	Bag volume	Bottom Outlet (free flow)	NRV in bag
84077181	2000 ml	Yes	Yes
84062182	3000 ml	Yes	Yes

Material specifications:

Contains latex or DEHP?	No
Bag incl. bottom outlet	PVC – DEHP free, Versapor RC membrane, PE, ABS
Packaging materials	PE, PET/PE composite film, Corrugated cardboard

Certifications:

CE-marking?	Yes (CE 2862)
Product classification	Class Is (EU)
Quality system certified to	ISO 13485
Global Medical Device Nomenclature (GMDN)	32072
EMDN Code	A06030302

Sterilization:

Sterile?	Yes
Sterilized by	Ethylene oxide
Shelf life	3 years

Storage and handling

Instructions for use included.



Packaging specifications:

REF No	Quantity (pcs)	Dimensions (mm)	Gross Weight (kg)	Volume (m³)
84077181	50	400x360x200	3,5	0,029
84062182	40	400x360x200	3,6	0,029

